


Thank you for purchasing the Rexall[®] Jumbo Flex Digital Thermometer. To correctly utilize all the advantages of this product, please familiarize yourself with these operating instructions before using the thermometer for the first time. Keep the instructions in a safe place for future reference.

FEATURES:

- Memory: The last reading is automatically displayed for 2 seconds when the unit is switched ON. Up to 9 sets of memory can be stored.
- High Temperature Indication: Rapid beeps and a red backlight alert the patient that he/she may have a temperature above 37.6°C / 99.7°F.
- °C / °F switchable
- Measurement time: 8 seconds

SAFETY INFORMATION:

- High, prolonged fever requires medical attention, especially for young children. Please contact your physician should this occur
- This thermometer contains small pieces that can pose a choking hazard. Do not allow children to take their own temperatures unattended.
- Keep battery away from children, harmful if swallowed.
- Do not open the thermometer except to replace the battery.
- To avoid damage to the probe, do not bend beyond a 45° angle.



HOW TO TAKE A TEMPERATURE:

1. Select Oral, Axillary (underarm) or Rectal as your measurement site based on each individual. Clean thermometer with soap and water prior to use.
2. Push the operating button located on the top of the thermometer. The thermometer will beep and the display will read: as illustrated in **Figure 1**. This is a function check indicating that the unit is working properly. Immediately after, "Last memory and the temperature will be displayed, then "----" with the °C/°F flashing will indicate the thermometer is ready to use, as illustrated in **Figure 2**.
NOTE: Wait for the "°C" to flash 3 times before taking a measurement. The thermometer will beep twice indicating the thermometer is starting to measure.
3. While taking a temperature, the °C/°F symbol flashes as long as the temperature increases.
4. The thermometer will beep three times, indicating that the peak temperature has been reached, right after the beep is heard the screen will light up green if the temperature is below 37.6°C / 99.7°F. If the temperature is 37.6°C / 99.7°F and above, rapid beeps will sound and the screen will light up red indicating a fever.
NOTE: The LCD will display "HI°C" if the measurement taken is below 32.0°C / 89.6°F. It will display "HI°C" if the measurement is above 42.9°C / 109.2°F.
5. To prolong the life of the battery, turn the thermometer off after each use by pressing the operating button briefly, otherwise the unit will turn off automatically in approximately 3 minutes.
NOTE: For best results, wait 3 minutes between temperature measurements.
6. To switch the thermometer from °C to °F, press and hold the operating button when the thermometer is turned off, once "----" is displayed with the °C or °F not flashing, press the operating button again to toggle between °C and °F.



Figure 1



Figure 2

WHAT IS "NORMAL" BODY TEMPERATURE?

"Normal" body temperature varies from person to person and changes according to time of day and body location. The following table shows the temperature ranges:

	Fever Temperature	Temperature Ranges
Rectal*	Over 38°C / 100.4°F	36.6°C – 38°C / 97.83°F – 100.4°F
Oral (mouth)*	Over 37.5°C / 99.5°F	35.5°C – 37.5°C / 95.9°F – 99.5°F
Axillary (underarm)*	Over 37.2°C / 99°F	36.5°C – 37.2°C / 97.7°F – 99°F

*Fever definitions according to the Canadian Paediatric Society

IMPORTANT: It is recommended to measure the body temperature when healthy to establish the specific "normal" temperature.

WARNING! A high temperature requires seeking immediate medical attention.

ORAL USE:

- Do not drink hot or cold fluids, exercise, smoke or perform other activities prior to a measurement that will raise or lower temperature readings when compared to your normal average temperature. The mouth should remain closed for up to 5 minutes before attempting a reading.
- Place the probe tip well under the tongue and wait until a beep is heard. Instruct children not to bite down on the probe.

AXILLARY (UNDERARM) USE:

- Wipe arm with a dry towel.
- Place the probe tip in the hollow of the armpit and bring arm down against the body and wait until 3 beeps are heard.

RECTAL USE:

- Commonly used for babies, young children or when it is difficult to take an oral or underarm temperature.
- Lubricate the thermometer with a water-soluble personal lubricant or petroleum jelly.
- Gently insert the probe tip no more than 0.5 inch or 1.3 cm into the rectum, and wait until a beep is heard.

DO NOT FORCE THE TIP INTO THE RECTUM IF RESISTANCE IS ENCOUNTERED.

MEMORY RECALL:

This thermometer recalls the last 9 temperatures. To obtain this information, press the power button to turn the thermometer on. When screen displays the last temperature taken, press the power button again to enter into memory recall. Keep pressing the button to view last 9 memories.

CLEANING AND DISINFECTING:

Before and after each use, clean the thermometer with warm soapy water and a soft cloth. Disinfect using a cotton ball moistened with alcohol. Never immerse in hot water or in harsh cleaning solutions.

CHANGING THE BATTERY:

When the **BAT** appears on the display during a measurement, the battery is exhausted and needs replacing. Remove the battery compartment located on the back of the thermometer, use a coin or flat head screwdriver to open. Discard the old battery and insert a new equivalent battery (CR2032, 3V). Ensure the battery is correctly positioned with the "+" facing up or towards the back of the unit. Replace battery compartment lid.

LIFETIME WARRANTY:

This product is warranted to be free from defects in material or workmanship for the life of the original owner. It will be repaired or replaced without charge upon receipt of the unit prepaid with \$5.00 to cover handling, packaging and return postage to the address below. Please include your full name, address and daytime phone number or email address.



Read the instruction manual carefully before using this device, especially the safety instructions, and keep instruction manual for future use.



Type BF Applied Part



Batteries and electronic devices must be disposed of in accordance with the local applicable regulations, not with domestic waste.

IP22: Protected against solid foreign objects of 12.5 mm diameter and greater. Protected against vertically falling water drops when the device is tilted up to 15°.

Quality Assured / Qualité Assurée
 Manufactured for / Fabriqué pour :
 Rexall Pharmacy Group Ltd.
 Mississauga, ON L4Z 1R9 ©2017
www.Rexall.ca
 by / par: Thermor Ltd.
 16975 Leslie Street, Newmarket, ON L3Y 9A1

BPA Free

