



Be. Better™/MC Blood Pressure Monitor Comparison Guide



Be. Better™/MC
Premium Plus Blood
Pressure Monitor

Name and Model Number	Automatic Wrist Blood Pressure Monitor RX803	Automatic Blood Pressure Monitor RX808	Premium Blood Pressure Monitor RX805	Premium Plus Blood Pressure Monitor RX809
Image				
Features				
Automatic Inflation	✓	✓	✓	✓
Time and Date	✓		✓	✓
Number of Users	1	1	2	2
Memory	90	30	60	120
Averages Last 3 Readings or Averages All Stored Readings	Averages Last 3 Readings	Averages Last 3 Readings	Averages Last 3 Readings	Averages All Stored Readings
Validated by Hypertension Canada		✓	✓	✓
Cuff Size Included	12.5 - 20.5 cm	Wide Range Cuff (M/L) 23 - 43 cm	23 - 36 cm	Wide Range Cuff (M/L) 23 - 43 cm
Batteries Included	✓	✓	✓	✓
Adapter Included		✓	✓	✓
Storage Bag Included	✓	✓	✓	✓
Clinically Validated		✓	✓	✓
Irregular Heartbeat Detection	✓	✓	✓	
Hypertension Indicator	✓	✓	✓	✓
Blood Pressure Analyzer Software				✓
Patented Pulse Arrhythmia Technology				✓
Warranty	Lifetime	Lifetime	Lifetime	Lifetime
UPC	77105840168 1	77105840476 7	77105840167 4	77105840475 0

† **European Society of Hypertension International Protocol / † Société européenne d'hypertension Protocole international**

The "International Protocol" also known as the EHS2010 protocol (European Hypertension Society) is an internationally recognized methodology used to determine the accuracy of a home use blood pressure monitor.

‡ This monitor meets Hypertension Canada's accuracy standard according to the International Protocol for automatic blood pressure devices. The protocol outlines specific requirements regarding the sample of subjects, test equipment and how the measurements are to be taken.

The studies are undertaken in clinical settings and the results are published in the Journal of Blood Pressure Monitoring and come under peer review. Devices that have met all the accuracy conditions are then considered "validated".

Helps detect an irregular or erratic heartbeat, medically referred to as an arrhythmia, can be intermittent and therefore difficult for doctors to detect. If the PAD icon appears frequently during measurements consult your doctor.



**Blood Pressure Hotline:
1-866-536-2289**

**Available exclusively
at Rexall**